

DKK-TOA CORPORATION

AUT-801 is a titrator that meets a variety of applications, such as simultaneous two-line titration.

Simultaneous titration of 2 series

By adding an optional burette to the standard configuration (single system), two different titration such as neutralization and oxidation-reduction titration can be performed simultaneously. In addition, it can be used for various purposes such as simultaneous measurement of pH and neutralization titration. In addition, two sets of multi-sample measurement devices (turntables) are connected, and can construct a two-line multi-sample simultaneous titration system.

Single System

A system without simultaneous based on standard configuration

Dual System



Dual-strain simultaneous titration system

In addition to the standard configuration, a titration burette (ABT-8), stirrer (ST-8), electrode, electrode stand, and electrode holder are required separately.

Supports simultaneous display of two lines



Simultaneous display and titration of system A and B are supported. Viewing data and reanalyzing is available from one system while measuring with another.

Versatile functions to support titration

Support functions including optimization of titration conditions

Titration time shortened with an operation similar to manual analysis (semi-automatic mode)





Compatible with various titrations lineup of high-performance electrodes

- ◆ The pH electrode uses a pH glass electrode (Strong-pH) that is hard to break (excluding some electrodes).
- We offer various types of electrodes suitable for various titration methods, such as ORP electrodes for redox titration, silver composite electrodes for salinity titration, electrodes for non-aqueous titration, and EC cells for electrical conductivity titration.

Compatible with various titrations

Equipped with 2 channels of pH/mV input as standard. Optional titration units (light intensity, polarization, electrical conductivity, potential difference) can be added to the 3rd and 4th channels.

Data Memory Function

At the time of single system, 600 data can be memorized in the main body.



Data Print Example

System [A] [Titration Resul Sample No.002 Date 20 Titration Time Operator Name Sample Size	1t])23/ 2/ 13 11:50 00:01:23 GUEST 1.0000
Mode No. Neutral.tit.1 Electrode Ch. Buret Ch. Wait Time Tit.Step Tit.Unit Valid Diff.1 Cont.P.1(Unit) Over Tit.Cnt.1 Fast Tit.	1 1 0 1 PH(ATC) 400.0 NoSet 4 Normal

The built-in printer prints in English, Korean, and Japanese making it easier to understand the titration conditions.

Equipped with validation support function

PH calibration history (up to 20 calibrations), burette capacity verification history creation function.



Potentiometric Titration Unit (FUT-8050) Up to 2 types connectable

(2 types of same unit are connectable)



System Example

Multi-sample Automatic Titration System AUT-801(single) + TTT-710

By setting the pretreated sample solution (beaker) on the turntable, titration of multiple samples is automatically performed. By using the flex function, it is possible to set the measurement conditions for each sample.

*The titration burette can also be connected to the previous model ABT-7.

Dual Multi-sample Automatic Titration System AUT-801(dual system) + TTT-710×2

Multi-sample titration is available for two different types of measurement item.

Also it can be utilized for upgrading the titration process function of the same samples.

*Titration burettes can also be connected to the previous model ABT-7.

Multi-sample Automatic Sampling Titration System AUT-801 + TTT-710 + AST-3210

By setting an appropriate amount of sample (beaker) on the turntable in advance, a series of processes such as sample weighing, dilution, reagent dispensing, titration, and washing can be automatically performed for multiple samples.





Auxiliary Equipments

Supports continuous automatic measurement of multiple samples

Turntable **TTT-710**

Excellent Maintenance

Electrode cleaning tank and electrode storage tank placed on the front of the device. Maintenance of pipes, electrodes, etc. has become easier.

Abundant electrode cleaning mode

Equipped with shower cleaning with pure water as standard. Chemical cleaning, bubbling cleaning, and air blowing are available as options.

Equipped with a stand for the analyzer (optional)

In order to save space including the analyzer, we have prepared a dedicated analyzer installation stand.



Specification

12 samples		200mL beaker, 300mLtall beaker (Recommend: AGC Techno Glass or HARIO glass beaker)		
Number of 18 samples	100mL tall beaker (Recommend: AGC Techno Glass or HARIO glass beaker)			
samples/ 36 samples response beaker 60 samples		30mL beaker, 50mLtall beaker(Recommend: AGC Techno Glass beaker)		
		20mL dedicated beaker (Nichiden Rika Glass Co., Ltd. H-20)		
	100 samples	20mL dedicated beaker (Nichiden Rika Glass Co., Ltd. H-20)		
St	irring method	Standard:Magnetic stirring method Optional:Propeller Stirring method		
Cleaning mode		 ① Pure water shower → Air blow ② Pure water bubbling → Pure water shower → Air blow ③ Chemical shower → Pure water shower → Air blow ④ Chemical bubbling → Pure water shower → Air blow ⑤ Chemical bubbling → Pure water bubbling → Pure water shower → Air blow ⑥ Chemical bubbling → chemical shower → air blow ⑦ Chemical shower → Air blow ⑧ Chemical bubbling → Pure water bubbling → Chemical shower → Air blow ⑦ Chemical bubbling → Pure water bubbling → Chemical shower → Air blow ⑧ Chemical bubbling → Pure water bubbling → Chemical shower → Air blow ⑧ Chemical bubbling → Pure water bubbling → Chemical shower → Air blow Note) •Air blow cleaning is available in the optional air pump box at time of usage. •Bubbling cleaning is available with optional air pump box and waste liquid valve at time of usage. 		
Cleaning time		0 to 9999 sec.		
End s	sample detection	Detecting by end detector pin or end table number setting		
Clean	ing Tank Material	PP		
Performance guarantee temperature		5 to 40°C		
	Power	AC100 to 240V 50/60Hz		
Powe	Standard	AC100V at time of usage:Max. Approx.60VA AC240V at time of usage:Max. Approx.90VA		
Supply Optional Co	Optional Connection	AC100V at time of usage:Max. Approx.100VA AC240V at time of usage:Max. Approx.130VA		
Dimensions		Approx.440(W)×409(H)×391(D)mm (When the table plate is not attached)		
		Approx.566(W)×409(H)×507(D)mm (When the table plate is attached)		
		Approx.566 (W)×534 (H)×507 (D)mm (During operation Max.dimension)		
Main Unit Weight		Approx.16kg (When the table plate is not attached) Approx.18kg (When the table plate is attached)		



Standard Accessories

Table (Select one from 12, 18, 36, and 60 samples)	End detecting pin
Table (Select one from 12, 10, 30, and 00 samples)	O ring
Electrode cartridge	Small hoffman pinch cock
(Numbers of samples/specify one type depending on the application)	Power code
Cleaning water tank (with 10L solution level sensor)	2P·3P transmitter adapter
Waste solution tank (with 10L solution level sensor)	Disposal beaker(200ml)
Cleaning water tube (3m)	Instruction Manual
Waste water hose (1.5m)	

Option

Product name	Code
For connecting turn table RS-232C cable (2m) (For connecting AUT-801)	7703820K
Air pump box (For air blow)	7400560U
Air pupm box (For air blow + bubbling)	7401640U
Waste liquid valve for TTT-710 (Pinch cock)	7401650U
Waste liquid valve for TTT-710 (Solenoid valve)	7401660U
Propeller stirring unit (For 12, and 18 sample table)	7401670U
Propeller stirring unit (For 36 samples)	7401680U
Propeller stirring unit (For 60, and 100 sample table)	7401690U
Installation table for analyzer (For AUT-801)	7401710U

Electrode Cartridge

Product name (Numbers of connectable electrodes / nozzles)	Code	Exterior		Product name (Numbers of connectable electrodes / nozzles)	Code	Exterior
1CH cartridge 1 (X) (X mounting hole:1, N mounting hole:1) (12, 18, 36, 60, 100 samples)	7505010K	Electrode Nozzle		2CH cartridge 2(X) (X mounting hole :2) (12, 18, 36samples)	7505030K	Electrode
1CH cartridge 2(X) (X mounting hole:1, N mounting hole:2) (12, 18, 36 samples)	7505020K	Electrode Nozzle		2CH cartridge3 (X mounting hole :2, N mounting hole :2) (12, 18samples)	6597940K	Electrode Nozzle
1CH cartridge3(X) (X mounting hole:1) (12, 18, 36, 60, 100samples)	7506840K	Electrode		Multi cartridge (J mounting hole :2, ION mounting hole :2, N mounting hole :4) (12samples)	6597980K	Nozzle
2CH cartridge1 (J mounting hole :2, N mounting hole :1) (12, 18, 36samples)	6597970K	Electrode Nozzle		Sample suction cartridge (N mounting hole :1) (12, 18, 36, 60, 100samples)	6597990K	Nozzle

*X, J, ION, N indicated on the mounting hole describes mounting electrodes and nozzels. X:for X/R series electrode J:for J/S series electrode ION:single-purpose type ion electrode

*Beakers for each table plate

N:nozzels

12 sample table plate :200mLbeaker, 300mLtall beaker
36 sample table plate :30mLbeaker, 50mLtall beaker
18 sample table plate :100mLtall beaker
60 sample table plate :20mL Dedicated beaker

● 60 sample table plate :20mL Dedicated beaker

6



Automation of measurement to titration to clean

Automatic sampling measuring equipment
AST-3210



Detector body

Sample lightweight method	Three-method cock switching measuring tube system
Amount of sample	Approx. 0.2 to 20mL (loop volume (measuring tube) fixed)
Light weight repeatability	Not more than 0.5% CV (with sample volume 10mL)
	Material : Hard glass
	Maximum capacity : Approx. 150mL
Measuring tank	Minimum requirement : Approx. 100mL
	Number of Electrodes Mounted : Up to 4 bottles
	Number of titration nozzles installed : Up to 4 bottles
Alarm display	Pure water empty, waste liquid, system failure
Waste liquid/pure water porin tank	20L polyethylene tank
Power supply	AC100V 50/60Hz
Power consumption	Approx. 100VA
Dimensions•Weight	Approx. 565 (W) × 635 (H) × 480 (D) mm Approx. 45kg

*Connection cables to the titrator main unit, electrode relay leads, etc. are required separately.

Optional Detection Unit

Photometric titration unit FUT-8010

(photometric sensor probe OPE-21A attached)

Used to detect color change by indicator method. It can be used for an appropriate amount of metal ions such as plating solution.

*Adapter is incorporated

into the main body, and

detector main body is

separately placed.



as ha

Detector body dimensions : Approx. 102 (W) \times 153 (H) \times 119 (D) nm Standard-Interference Filter : 530nm, 630nm

Adapter



Measurement method : Constant current voltage method or constant voltage current method Applied voltage : 0 to ±1000mV (arbitrarily set) Applied current : 0 to ±99.9µA (arbitrary setting)



CONDUCTIVITY CEL

The unit is integrated into the main body.

Measuring range : 0 to 200.05/m (depending on the cell used) Manual range for the titration mode (The range is set by the main unit) Temperature : 0 to 100.0°C

Small titration unit

Small titration-kit ASSY **7075600K**

We have prepared a small titration unit combined with a pH composite electrode for trace quantities, a silver composite electrode, etc.

The amount of titrant waste solution can be reduced.

(Some titrations may not be available.)

Min. volume : Approx. 10mL Max. volume : Approx. 20mL

Configuration

Adapter, stirrer (2 pcs), Titration nozzle Measuring cell (50 pcs.), O-ring (P5) (2 pcs.), O-ring (P12) (2 pcs.)



Electrodes

Uses (Representative titration examples)	ltem name	Model/ Code No.	Note
	pH compound electrode (For general use)	GST-5841C	AUT-801 standard-attached electrode
Neutralization	pH Combined Electrode (For Trace)	GST-5845C	
titration	pH combined electrode (Double junction)	ELP-062	Liquid junction replaceable Useful when spillage of KCL solutions is problematic
	ORP compound electrode (For general use)		
Oxidation reduction titration	ORP composite-electrode (Double junction)	ELM-016	Liquid junction replaceable Useful when spillage of KCL solutions is problematic
Salinity titration	Silver composite electrode (Double junction)	ELX-006	Liquid junction replaceable
Nonaqueous	Nonagueous Glass electrodes (For general use)		Electrode adapter (0JD00001) for connecting the main unit
titration	Comparison electrode (Double junction)	HS-305DS	An electrode-holder (S-HLD-S) is required separately
Photometric titration	Light sensor probe	OPE-21A	Photometric titration unit standard attached probe
Polarizable titration	Two platinum electrodes	HPD-303	Polarized titration unit standard attached electrode
Titration electrical conductivity	Conductivity cell	CT-58101B	Electrode with electrical conductivity titration unit standard

Use of pH glass electrodes that are resistant to cracking

We have used Strong & Float-pH electrodes that combine our sensor technologies. (excluding some electrodes)

Strong UP the strength of the tip (more than 10 times that of us) Float Built-in float that can be checked for inner liquid density at a glance



8

Parts and standard solutions

Main unit related parts Standard	Accessories: Additional a	rrangements are required when constructing the standard accessory dual system.
ltem name	Model/Code No.	Note
Stirrer	ST-8	
Electrode stand * Standard Accessories	7702590K	With support and stopper
Electrode Holder Standard Accessories	7430850K	
Electrode attachment (G) **Standard Accessories	0IB00004	For electrodes of X series etc.
Electrode attachment (J) **Standard Accessories	0IB00005	For electrodes of J series etc.
Electrode attachment (N) Standard Accessories	0IB00008	For nozzle, etc.
Electrode attachment (ION)	0IB00006	Single capacity ion electrode
Electrode adapter	0JD00001	Single-function electrodes can be connected
Electrode holder (R series)	0IB00001	Required when attaching 4 to 5 electrode attachments
Electrode relay lead (1m)	0GB00001	
Electrode relay lead (3m)	0GB00002	
Minor titration kit	7075600K	Configuration: Adapter, Stirrer (2 pcs.), Titration Nozzle, Measuring Cell (50 pcs.), O-ring (P5) (2 pcs.), O-ring (P12) (2 pcs.)
Turntable Connecting Cable (2m)	7703820K	For connecting TTT-510/TTT-710
AST connecting cable (2m)	7703830K	For connecting AST-3210
RS-232C connecting cable (2m)	0GC00002	For PC connection (PC connector standard: D-Sub9P)
USB Communication Cable (2m)	7473100K	For connecting a PC (cable length 2m)
Comparative electrode-internal-liquid RE-1 100mL	143F230	Internal solution of monofunctional comparative electrode (HS-305DS, etc.)
Comparative electrode-internal-liquid RE-2 100mL	143F238	Outer cylinder liquid such as an electrode for salinity titration
Printer paper (Volume 5)	PAP-HCS	AUT-801 Standards Attached/Thermosensitive Paper for Built-in Printer
OA tap out WCH2436H	102DE26	If the number of outlets is insufficient, purchasing is required Wire length 3m, 6 ports

Bullet-related parts

Item name	Model/Code No.	Note
Syringe unit base (for 1 to 20mL syringes)	7703030K	Standard accessories for ABT-8 Syringes, nozzles, tubes, and reagent bottles are sold separately %If a syringe other than "20mL syringe module brown" is used, a syringe holding ring is required separately
Syringe unit base (for 50mL syringes)	7708720K	Dedicated for ABT-8 : Syringes, nozzles, tubes, and reagent bottles are sold separately
Syringe holding ring (syringe (brown) 1 to 20mL)	X8774900	
Syringe holding ring (for syringe (brown) 50mL)	X5009600	
Cap nut (for 3-way cock) 5pcs	AUT-061P	
Contains 5 bag nuts (for syringe)	AUT-062P	
Tapered piece (10 pcs.)	AUT-063P	
20mL Syringe Unit (brown)	7702900K	Standard accessories 20mL syringes (for modules), syringe covers (brown), syringe caps for 20mL, bag nuts, and one tapered piece



Bullet-related parts

Item name	Model/Code No.	Note
For 20mL Syringe Module Unit	7702910K	20mL syringe-unit (brown) components
Syringe cover brown	76808800	20mL syringe-unit (brown) components
20mL syringe cap	76274500	20mL syringe-unit (brown) components %O-ring for 20mL syringes, syringe packing, bag nut, and tapered piece are sold separately
O-ring for 20mL syringes	115A089	20mL syringe-unit (brown) components
20mL syringe packing	115J074	20mL syringe-unit (brown) components
Retaining washer	77078700	1 bottle For retaining pipes
50mL syringes brown	P000010	With one cap nut and one taper piece For ABT-8 use, a separate syringe-unit base (50mL) is required
20mL syringes brown	AUT-045P	With one cap nut and one taper piece For ABT-8 use, a separate syringe-holding ring (X8774900) is required
10mL syringes brown	AUT-046P	With one cap nut and one taper piece For ABT-8 use, a separate syringe-holding ring (X8774900) is required
5mL syringes brown	AUT-047P	With one cap nut and one taper piece For ABT-8 use, a separate syringe-holding ring (X8774900) is required
1mL syringes brown	AUT-048P	With one cap nut and one taper piece For ABT-8 use, a separate syringe-holding ring (X8774900) is required
Defoaming nozzle (for 1 to 20mL syringes)	P000070	Standard accessories Brown with bag nut
Defoaming nozzle (for 50mL syringes)	P000071	Brown with bag nut
Reagent bottle (with tube joint)	7075630K	ABT-8 reference product attached 1000mL semi-transparent brown poly bottle
Carbon dioxide absorption tube (with tube)	7075640K	ABT-8 reference product attached
Teflon tube black (for 5-50mL syringes) 2m	AUT-022P	ABT-8 reference product attached For 5 to 50mL syringes (ϕ 2 x ϕ 3)
Teflon tube black (for 1mL) 2m	AUT-024P	For 1mL syringes only (ϕ 1 x ϕ 3)
Hardware for syringes	AUT-066P	ABT-8 reference product attached
Cleaning metal fitting for nozzle	AUT-067P	ABT-8 reference product attached
Nozzle purge tube	P000108	ABT-8 reference product attached
Adapter	70774700	Minor titration kit components
Stirrer bar	107D101	Minor titration kit components
Titration nozzle (for trace amounts)	AST-P008	Minor titration kit components
H-20(50 pcs. / 1 box)	136C591	Minor titration kit components
O-ring (P12)	115A620	Minor titration kit components
O-ring (P5)	115A628	Minor titration kit components

Titrator body AUT-801

Specifications			
Display	7inch color graphics LCD		
Display language	Japanese / English / Korean		
Display content	Titration curves, titrations, pH or mV, liquid temperature, and titration times concurrently displayed For dual-system: System A/ system B simultaneous display/switch display is enabled		
Operation key/ operation method	Flat key/function key, numeric keypad interactive method Numeric keypad toggle input supported		
JIS Model (pH)	JIS type I		
Law Type Approval No. (pH)	Application pending		
	Inflection point detection titration	Stat titration	
	Setpoint detection	Pre-titration	
Titration /	(titration unit/derivative value)	pKa determination	
method	Total titration	pH/mV determination	
	(inflection point detection/set point detection)	pH control	
	Intersection detection titration	Electrical conductivity measurement	
Titration stage number	Up to 5 stages (intersection d	etection up to 2 stages)	
Titration mode	Standard mode: 20/System User mode: 40/system Link Mode (For Sequence Tit	ration): 20/System	
Titration control system	Single System (System A) Switching dual system (System A and System B) sets Dual-system simultaneous two-series titration/measurement enabled		
Measurement items/range	pH:pH 0.00 to pH 14.00 reference 2ch only mV:0.0 to ± 2000.0 mV Temperature :0.0 to 100.0°C [When the option is connected below] Galvanostatic method: 0-2.000 μ A 0-20.00 μ A 0 to 200.0 μ A 0 to 1000 μ A (depending on the range switching) Constant current/voltage method: 0.0 to ± 2000.0 mV Electrical conductivity: depending on the cell used Display range 0 to 200.0 μ S/m (0 to 2.000 μ S/cm) 0 to 2.000 mS/m(0 to 20.00 μ S/cm) 0 to 20.00 mS/m(0 to 20.00 μ S/cm) 0 to 20.00 mS/m(0 to 20.00 mS/cm) 0 to 20.00 mS/m (0 to 2.000 mS/cm) 0 to 20.00 mS/m (0 to 2.000 mS/cm) 0 to 20.00 mS/m (0 to 2.000 mS/cm) 0 to 20.00 m (0 to 2.000 mS/cm) 0 to 200.0 m (0 to 2.000 m m m m m m m m m m m m m m m m m m		
Number of electrode inputs	Max. 4ch Potentiometric (pH/mV)2ch equipped as reference) 3ch, 4ch is added by the optional detecting unit		
Number of titration burette connections	Up to 10 units can be linked		
pH Calibration	Automatic 5-point calibration (manual calibration f	or optional standard solution calibration)	
pH reference solution selection	JIS Standard Solution / US Stan Optional Standard Solution (Ma	dard Solution / Type 2 / x. 5-point Calibration)	
Standard solution for pH optional	Input 2 arbitrary standard so	lution table	
pH Temperature compensation range	ATC (automatic temp. compensation): 0.0 to 100.0°C MTC (manual temp. compensation): 0.0 to 100.0°C		

Temperature ele conductivity corr	The temperature compensation range	ATC (automatic temp. compensation): 0.0 to 100.0°C MTC (manual temp. compensation): 0.0 to 100.0°C No Temp. Compensation (ATC OFF)		
ctrica	Reference temperature setting	0 to 100.0°C		
lation	Temperature coefficient (linear)	0 to 10.0%/°C		
Data (re-a	a memory analyzeable)	Up to 600 data per series $300 \times 2 = 600$ data for multitasking in 2 series (dual system)		
Print	ter	Built-in line thermal printer		
Validation Support Functions		pH Calibration History: Up to 20 Calibration minutes Voluntary Inspection History: Up to 6 (Voluntary Inspection with Checker Input) Periodic Inspection History: Up to 10 (Periodic Inspection With Standard Reagents) Bullet Capacity Assay History: Up to 6 (The burette has information) Instrument and electrode control deadline alarm (date control) Reagent Replacement Alarm (Date Control) Syringe Replacement Alarm (Select either Date Control/Stroke Count Control) Reagent remaining amount warning (set reagent amount beforehand)		
Con func con	nmunication tion/external trol input/output	RS-232C×5ch 1.Titration burette (up to 10 units can be linked) 2.PC (D-sub) 3.Turntable (TTT-710) 4.Any one of Automatic sampling measuring device (AST-3210), Turntable(TTT-710) and electronic balance 5.External printer (for plain paper printing) (EPS-P30) 1 USB (for memory) 1 USB peripheral		
Other functions		Burette connecting function		
Perfo temp	ormance guaranteed perature and humidity	5 to 35°C 20 to 85% (non-condensing)		
Pow	ver supply	AC100 to 240V 50/60Hz (dedicated AC adapter)		
Pow	ver consumption	Max. 45VA(AC100V) Max. 60VA(AC240V) (Adapter DC Out DC24V 33VA)		
Dim	iensions	Approx. 150(W) × 215(H) × 385(D) mm		
We	ight	Approx. 3.2kg		

AUT-801 reference product attached

Item name	Model/Code No.
Stirrer (with one stirrer)	ST-8
Stirrer-burette connection cable	7692410K
pH electrode (Strong-pH compound electrode)**	GST-5841C
pH6.86 reference solution 500mL*	143F192
pH4.01 reference solution 500mL*	143F191
Solution within the reference electrode (50mL)	-
Electrode stand (with support, stopper)	7702590K
Electrode holde	7430850K
Electrode Attachment (Type G)	0IB00004
Electrode Attachment (Type J)	0IB00005
Electrode Attachment (Type N)	0IB00008
Printer paper (Volume 5)	PAP-HCS
Power cord	118C252
AC adapter	134L070
2P-3P converter adapter	118C504
Ground wire (2m)	XL600697
Ring instructions manual	_

*When no electrode or standard solution is specified, it is not attached.



Titratable burette ABT-8

Specifications

Display	Syringe size indicators: LED indication Error indicators: LED indication	
Syringe size setting	Flat key	
Syringe	Syringe inner cylinder precision polishing type 20mL clear glass syringe with brown syringe cover Optional syringe (brown glass syringe) 1 mL, 5 mL, 10 mL, 20 mL, 50 mL	
Material of wet part	PTFE, PCTEF, hard glasses, fluorine rubber	
Tube used	Tube to be used ϕ 2 (inside) x ϕ 3 (outside) Black Teflon tube (ϕ 1 (inside) x ϕ 3 (outside) black Teflon tube when using 1mL syringe)	
Aspiration rate	Full stroke Approx. 20 sec.	
Discharge speed	50mL syringes: 2 to 150 mL/ min 20mL syringes: 0.6 to 60.0 mL/ minutes 10mL syringes: 0.3 to 30.0 mL/ 5mL syringes: 0.2 to 15.0 mL/ minutes 1mL syringes: 0.03 to 3.00 mL/ minutes	
Minimum output	50mL Syringes: 0.0025 mL 20mL syringes: 0.001 mL 10mL Syringes: 0.0005 mL 5mL Syringes: 0.00025 mL 1mL Syringes: 0.00005 mL	
Bullet accuracy	In 20mL syringes Total volume error: ±0.1 % Dispensing repeatability: ± 0.01 mL	
Validation Support Functions	Burette Volume Test History: Up to 6	
Power supply output for stirrer	DC5V for dedicated stirrers	
Performance guaranteed temperature and humidity	5 to 35°C 20 to 85% (non-condensing)	

Power supply	AC100 to 240V 50/60Hz (dedicated AC adapter)
Power consumption	Max. 35VA(AC100V) Max. 45VA(AC240V) (Adapter DC Out DC24V 28VA)
Dimensions	Approx. 113 (W) \times 396 (H) \times 348 (D) mm
Weight	Approx. 5.2kg

◆ABT-8 reference product attached

Item name	Code No.
Syringe unit base (for 1 to 20mL) Including test tubes for nozzle standing	7703030K
Reagent bottle (1000mL semi-transparent brown poly-bottle)	7075630K
Carbon dioxide absorption tube (with tube)	7075640K
20mL Syringe Unit (brown)	7702900K
Syringe mounting brackets	AUT-066P
Defoaming nozzle (for 1 to 20mL)	P000070
Cleaning metal fitting for nozzle	AUT-067P
Purge tube for nozzle	P000108
O-ring for syringe	115A089
Teflon tube black (2m)	AUT-022P
Burette connecting cable (1.8 m)	118B129
Retaining washer	77078700
Power cord ※	118C252
AC adapter *	134L070
2P-3P converting adapter ※	118C504
Earth wire (2m) ※	XL600697
Ring instructions manual	-

*No service outlets. One power supply is required for each burette body. Buy optional power taps (6 units) (Item Code: 102DE26) as needed.

Stirrer ST-8

Compliant beaker	Up to 200mL beaker	
Connector terminal	Power supply for upper agitation unit	
Designed	DC5V (dedicated cable/dedicated AC adapter	
Power supply	AC100 to 240V 50/60Hz optional)	

Standard accessories

Item name	Code No.
Stirring bar (Φ8×25 mm)	107D039

Power consumption	Max. 5VA(AC100V) Max. 9VA(AC240V) (Adapter DC Out DC5V 5VA)
Dimensions • Weight	Approx. 110(W)x73(H)x135(D) and Approx. 0.8kg

**Dedicated AC adapter instead of stirrer bullet connecting cable if purchased separately (7430880K) And the power cord (118C229) is included. If the dedicated AC adapter is used, there is no service outlet.



Analysis Solutions in Wide Range of Applications



Food Processing

Type of titration	Used electrodes (unit)	Titration content
Acid-base titration	pH multiple electrode GST-5841C	Fruit beverages by set point titration Acidity analysis of yogurt, etc.
		Continuous analysis of acidity and formol nitrogen in fruit beverages
		Analyze the emulsion ($C_3H_6O_3$) concentration of the noodle processing solution
	pH multiple electrode GST-5823 S	Total acidity and amino acid content analysis of sake
Oxidation-reduction titration	ORP compound-electrode PST-5821C	Vitamin C concentration analysis in juices
		Peroxide value (POV) analyses of edible oils
Precipitation titration	Silver combined electrode ELX-006	Seasonings (soy sauce, sauce, vinegar, etc.) Salinity analysis of dried sauce
Chelatometric titration	Photometric titration unitFUT-7010 (filter 630nm)	Calcium concentration analysis in cornflakes
	Calcium ion-electrode CA-135B Reference electrode HS-305DS	Calcium concentration analysis in nutritional fortification agents
Nonaguagus titration	Glass electrode HGS-2005 Reference electrode HS-305DS	Oxidation analysis of edible oils
Nonaqueous titration		Analysis of Amino Acids by Perchlorate Titration Method
Titration electrical conductivity	Conductivity titratable unit FUT-8040	Sulfate ion (SO $_{4^{2-}}$) concentration analysis
Acid-base titration Precipitation titration	KCL Supply-type pH combined electrode ELP-062	Continuous analysis of acidity and salinity of
	Silver combined electrode ELX-006	vinegar, mayonnaise, etc.



Chemical and Analytical

Type of titration	Used electrodes (unit)	Titration content
Acid-base titration		Phosphate (H ₃ PO ₄) concentration
		Phosphate ester concentration analysis
	pH multiple electrode GST-5841C	With barium(Ba(OH) ₂ hydroxide.Fractional and appropriate amounts of barium carbonate (BaCO)
		By the sodium sulfite method Formalin (HCHO) concentrations analyses
		Hydrazine (N ₂ H ₄) concentration-analysis
		Analysis Iodine (I ₂) Concentration
		Hydrogen peroxide (H ₂ O ₂) concentration analysis
		Sodium Sulfite (Na ₂ SO ₃) Concentration Analysis
		Purity Analysis of Sodium disulfite (Na ₂ S ₂ O ₅)
Oxidation-reduction	OPP compound-electrode PST-5821C	Purity of ammonium persulfate ((NH_4) ₂ S ₂ O ₈)
titration	OKF Compound-electrode F31-3621C	Potassium Ferricyanide (K ₃ (Fe(CN) ₆]) Concentration Analyses
		In hydrogen sulfide absorption caustic soda
		Analyses of Flowing Soda (NaSH)
		Sulfur-ion (S ²⁻) densitometry
		Hydrazine (N ₂ H ₄) concentration-analysis
	Chloride ion-electrode CL-125B Reference electrode HS-305DS	Salinity analysis in raw concrete
Precipitation titration		Chloride ion (Cl ⁻) concentration (JIS R 5202) In Portland cement
	Silver combined electrode ELX-006	Analysiss of low chlorine (Cl ⁻) levels
		Sulfur-ion (S ²⁻) densitometry
Chelatometric titration	Photometric titration unit FUT-8010 (filter 530nm)	In cobalt(CoCl2) chloride Cobalt ion concentration analysis
	Silver ion electrode CU-125 Reference electrode HS-305DS	In basic copper base (CuCl ₂ · Cu(OH) ₂) Copper ion concentration analysis
Nonaqueous titration	Class electrode HCS 2005	Fractional determination of fluorine (HF) and nitric acid (HNO ₃)
	Reference electrode HS-2005	Analysis of epoxy equivalents
		Saponification number analysis
Titration electrical conductivity	Conductivity titratable unit FUT-8040	Sodium hydroxide (NaOH)Fractional Determination of Triethylamine ((C2H5)2N)



Plating

Type of titration	Used electrodes (unit)	Titration content
Acid-base titration	pH multiple electrode GST-5841C	Analysis of Borate (H_3BO_3) Concentration in Nickel Plating Solution
		In permanganate solution NaOH (NaOH) concentration analyses
		Separate determination of sodium hydroxide (NaOH) and sodium carbonate (Na $_2$ CO $_3$) in cadmium cyanide plating solution
		Chromate (CrO ₃) concentration in the plating solution
		Analysis of Sn^{2+} content in solder plated solution
Oxidation-reduction titration	Platinum-combined electrode PST-5821C	Analyzing Sodium Hypophosphite (NaH ₂ PO ₂) In nickel sulfamate plating solution
		Rochelle salt (KNaC_4H_4O_6) concentration-analysis In the bronze plating solution
Precipitation titration	Combined silver-electrode ELX-006	Nickel-chloride (NiCl02) concentration In the nickel plating solution
		In the bronze plating solution Sodium cyanide (NaCN) concentration
	Photometric titration unit FUT-8010 (filter 530nm)	Nickel (Ni) concentration-analysis in nickel-plating solution
		Lead(Pb) levels in solder plating solution
		Zinc Oxide (ZnO) Concentration in Zinc Plating Solution
Chelatometric titration	Photometric titration unit FUT-8010 (filter 630nm)	In cadmium blue plating solution Cadmium(Cd) level analysis
	Copper-ion-electrode CU-125 Reference electrode HS-305DS	Copper (Cu) concentration-analysis in bronze plating solution
		Nickel (Ni) concentration-analysis in nickel-plating solution (Reverse Liquid Constant Method with Copper Ion Electrode)
Acid-base titration Chelatometric titration Precipitation titration	Photometric titration unit FUT-8010 (filter 530nm) pH Electrode ELP-062 Combined silver-electrode ELX-006	Borate (H ₃ BO ₃) in nickel plated liquid, All nickel (Ni) and nickel chrolide (NiCl ₂), Sequential analyses of nickel sulphate (NiSO ₄)

Type of titration	Used electrodes (unit)	Titration content
	pH multiple electrode GST-5841C	Acid concentration analysis in the etchant
		Fractional Determination of Total Acid and Iron in Stainless Steel Treatment Solution
		Determination of Sulfuric Acid (H ₂ SO ₄) and Copper (Cu) Concentrations in Print Board Etching Solutions
Acid-base titration		Analysis of Free Hydrochloric Acid Concentration in Mask Treatment Solution
		Acid-concentration, aluminium(Al)-concentration analyses in chemical conversion agents
		Ammonium fluoride (NH ₄ F) concentration Hydrofluoric (HF) in the etchant
		Densitometry of oxalic acid ((COOH) ₂) in an etchant
		Copper-ion (Cu ²⁺) densitometry (JIS M 8121)
		Chromium ion (Cr6+) densitometry
	Platinum combined electrode PST 5821C	Cobalt-ion (Co ²⁺) densitometry
Oxidation-reduction	Flatinum-combined electrode F3T-362TC	Analyzing Titanium(Ti ³⁺ , Ti ⁴⁺) Density in Etching Solution
titration		Analysis of Fe-ion (Fe ²⁺) Concentration in Pickling Solution
		Sodium Thiosulphate (Na $_{\rm 2}S_{\rm 2}O_{\rm 3}$) Concentration Analysiss In the desulfurization solution
	Platinum-combined electrode ELM-016	Analyzing the Concentration of ${\rm H_2O_2}$ in Chemical Polish Liquid
		Densitometry of oxalic acid $((COOH)_2)$ in an etchant
	Combined silver-electrode ELX-006	Chloride ion (Cl ⁻) concentration-analysis in fluxes
Precipitation titration		Thiourea ((NH ₂) ₂ CS) Concentration Analysis In the continuous surface treatment solution
		Concentration of thiocyanate(SCN ⁻) in desulfurized liquor
		Analyzing Silver Ion (Ag ⁺) Concentrations
		Nickel(Ni) concentrations in etchants
	Photometric titration unit FUT-8010(filter 530nm)	Purity-analysis of zinc-oxide (ZnO) (JIS K 1410)
Chelatometric titration		Lead-ion (Pb ²⁺) Density Analysiss
	Photometric titration unit FUT-8010(filter 630nm)	Analysiss of high-purity strontium chloride (SrCl $_2$) concentrations
	Copper-ion-electrode CU-125 Reference electrode HS-305DS	Copper(Cu) Concentration in Chemical Polishing Liquids
	Glass electrode HGS 2005	Acid value analysis of the flux
Nonaqueous titration	Reference electrode HS-2005	Hydrochloric acid (HCl) in the surface-treated solution, Fractional determination of hydrofluoric acid (HF)



Electricity, Steel and Metals Т

Applications



Environment







Pharmaceuticals, cosmetics, and Perfume

Type of titration	Used electrodes (unit)	Titration content	
Acid-base titration	pH multiple electrode GST-5841C	Alkalinity analysis of the upper and sewage	
		Acidity analysis of upper and sewage	
Oxidation-reduction titration	ORP compound-electrode PST-5821C	Residual Chlorine (Cl ₂) Level Analyses of Upper Water	
		Analyzing Potassium Permanganate (KMnO4) Consumption	
		Dissolved-oxygen (DO) concentration analyses of seawater	
Precipitation titration	Photometric titration unit FUT-8010	Low concentration sulfuric acid (SO4 ²⁻)	
	(filter 630nm)	concentration analysis (JIS K 0103)	
	Combined silver-electrode ELX-006	Chlorine-ion (Cl ⁻) analyses of surface water	
Chelatometric titration	Photometric titration unit FUT-8010 (filter 630nm)	Total hardness analysis of clean water	
		Calcium(Ca) in tap water,	
		Fractional determination of magnesium(Mg)	

Type of titration	Used electrodes (unit)	Titration content
Acid-base titration	pH multiple electrode GST-5841C	Control of Sodium Hydroxide (NaOH), Amines, and Potassium Carbonate (K_2CO_3) Levels in Desulfurized Effluent
Precipitation titration	Combined silver-electrode ELX-006	Control of H ₂ S levels in desulfurized effluent
		Chlorine-ion (Cl ⁻) analyses of desulfurized effluent
Jonaqueous titration	Glass electrode HGS-2005 Reference electrode HS-305DS	Acid value analyses of engine oil (JIS K 2501)
		Base Value Analysis (JIS K 2501) of Engine Oil
		Acid Number Analysis (JIS K 2101) of Electric Insulating Oil
		Analyzing Vinyl Acetate (C4H6O2)
		Analysis of the carbonyl value of fats and oils
	Polarizable titratable unit FUT-8030	Bromine Value Analysis (JIS K 2605) of Petroleum Products

Type of titration	Used electrodes (unit)	Titration content
Acid-base titration	pH multiple electrode GST-5841C	Sodium bicarbonate (NaHCO_3) concentrations in stomach drugs
Oxidation-reduction titration	ORP compound-electrode ELM-016	Purity-analysis of sodium sulfide (Na ₂ S)
Precipitation titration	ORP compound-electrode ELM-016	Analyzing Benzethonium Chloride(C_27H_42ClNO2) With sodium tetraphenylborate
Chelatometric titration	Photometric titration unit FUT-8010 (filter 530nm)	Analyzing Aluminium Oxide (Al2O3) Concentrations in Stomach Drugs
	Photometric titration unit FUT-8010 (filter 630nm)	Magnesium-oxide (MgO) concentration-analysis of stomach chemicals
		Analysiss of Precipitated Calcium Carbonate (CaCO ₃) Concentrations in Gastric Medicines
Nonaqueous titration	Glass electrode HGS-2005 Reference electrode HS-305DS	Ester value analysis of cosmetic perfumes
		Acid value analysis of cosmetic perfumes
		Cyclohexidin Gluconate (C ₂₂ H ₃₀ Cl ₂ N ₁₀) Concentration analysis
	Polarizable titratable unit FUT-8030	Ethyl Aminobenzoate(NH ₂ C ₆ H ₄ COOC ₂ H ₅) Concentration analysis
Titration electrical conductivity	Conductivity titratable unit FUT-8040	Sulphate (SO $_4^{2-}$) density analyses



Overseas Sales Division: DKK-TOA Corporation 29-10, 1-Chome, Takadanobaba, Shinjuku-ku, Tokyo 169-8648 Japan Tel : +81-3-3202-0225 Fax : +81-3-3202-5685 E-mail : intsales@dkktoa.com



Please read the operation manual carefully before using producuts.

https://www.toadkk.com/english