

# IA-300

## Ion analyzer

Easy simultaneous  
measurement of multiple ions

- ◆ Six cations or seven anions can be measured
- ◆ Anion are available for both suppressor and non-suppressor systems.



# Ion analysis is easier and more accessible!

The ion analyzer is a measuring instrument using ion chromatography as a measuring principle. Ion chromatography is a reliable method used in many official methods such as JIS. It is used not only in environmental measurement, but also in a variety of fields, including raw materials, production lines, quality inspections, and wastewater measurement. It is an ion analyzer to be used in more familiar areas.

## Ion analyzer

One packaging of ion chromatography functions (pump, sample injection, constant temperature, detection, data processing, printer, etc.) and injection of sample into the instrument can be used to determine the concentration of the purpose ion.


## For the measurement of

**Anion**  $F^-$ ,  $Cl^-$ ,  $NO_2^-$ ,  $Br^-$ ,  $NO_3^-$ ,  $PO_4^{3-}$ ,  $SO_4^{2-}$  simultaneously

**Cation**  $Li^+$ ,  $Na^+$ ,  $NH_4^+$ ,  $K^+$ ,  $Mg^{2+}$ ,  $Ca^{2+}$  simultaneously

(Note: simultaneous measurement of anions and cations is not possible.)

## Configuration of the device



# Features

## Set measurement conditions with mode selection

Six cations or seven anions can be measured using ion chromatography.

Measurement mode	Measured ion
1, 2 Divalent cation simultaneous measurement mode.	Lithium ion, sodium ion, ammonium ion, Potassium, magnesium, and calcium ions
Anion measurement mode	Fluoride ion, chloride ion, nitrite ion, Bromide, nitrate, phosphate, and sulfuric acid

## Automatic processing for analysis and calculation

Analysis and calculation of measurement data are performed automatically by the instrument itself.

Compared with ion chromatographs, which process data PC, the processing time can be drastically reduced.

Can be converted to hardness and nitrogen content.

(For details, refer to the measurement example on page 5.)

## A suppressor method was also adopted for anion measurement.

In addition to the conventional measurement conditions, an anion measurement mode using an ion exchange membrane suppressor has been added.

The measurement accuracy of phosphate ions is improved, and the measurement range of each ion is widened.

(For details, refer to the specifications on page 9.)

## Various reagents are available (sold separately).

Since various reagents are available for measurement, they can be used immediately without any special equipment such as pure water or balance.

## Reduced labor by automatic measurement of multiple samples


The autosampler ICA-700AS, sold separately, is available for automatic measurement of multiple samples.

(See page 6 for details.)


## Easy maintenance

All column changes and piping can be operated on the front.

It also provides on-screen, easy-to-understand support for maintenance and troubleshooting instructions.



Maintenance mode display example



Pressure check mode display example




# Various convenient functions


## Simple procedure for introducing the eluent

Small amount of leftover below 3mL even with in-line degasser.

The cock is loosened counterclockwise, and the attached syringe is used for priming.



Eluent displacement is also smooth.




Operate without disconnecting the piping

## Adoption of large isothermal parts

Columns, guard columns, suppressors, and detection cells were all housed in the thermostat.

Improved stability and reliability of measurement



```
Stand-by 14.23
Measure mode Anions2
Column PC1-205
PC1-205G
Eluent For PC1-205
Flow rate 1.1mL/min
Inj. vol 20uL
Std. sol IA-AS1

Start:ST/SP Data:DATA
Setup:FUNC
```

After starting the equipment, specific equipment status such as pressure and electric conductivity value is displayed.


## Timer function assuming various handling

You can set the time to start or stop the instrument for each day of the week.

```
Timer setup 4/25 Tue. 14.12
Wakeup Time Stop Time
Mon. ON 08:30 OFF 17:00
Tue. OFF 08:30 OFF 17:00
Wed. OFF 08:00 OFF 17:00
Thu. OFF 08:00 OFF 17:00
Fri. OFF 08:00 OFF 17:00
Sat. OFF 08:00 OFF 17:00
Sun. OFF 08:00 OFF 17:00

[Please choose ON/OFF/Weeks]
Change: ↑ ↓ Move: ← → Return:ESC
Enter:ENTER Copy:CAL End:FUNC
```

Measure immediately without waiting.





# Display and Printing Examples


Quantitative values and chromatograms can be displayed and printed.

## Anion calibration data

Example of calibration data display


Peak	Items	Values	Equivalents
1	P04	4.01 mg/L	
2	F	0.491 mg/L	
3	Cl	0.985 mg/L	
4	NO2	1.00 mg/L	
5	Br	0.992 mg/L	
6	NO3	0.983 mg/L	
7	S04	1.97 mg/L	

Chromatogram



### Printing example

This is an example of printing a calibration result. Calibration results are printed unconditionally regardless of the printer settings.




## Cation calibration data

Example of calibration data display


Peak	Items	Values	Equivalents
1	Li	2.01 mg/L	
2	Na	9.91 mg/L	
3	NH4	9.99 mg/L	
4	K	20.0 mg/L	
5	Mg	9.93 mg/L	
6	Ca	20.0 mg/L	

Chromatogram



### Printing example

This is an example of printing sample measurement results.




Printer error

No Paper:  
Use of the printer is stopped.  
Please settle the paper jam of printer  
and press the ENTER key.

Check:ENTER

Show chart paper refill  
message when chart  
paper is gone



Chromatograms can be  
printed on the specified  
display scale.


Printing to match the  
low density

# Example of measurement


## ●Example of conversion printing of nitrogen, phosphorus, sulfur, and hardness

### ■ Tap water

Anion		
Factor		
1 F	0.099	mg/L
2 Cl	16.0	mg/L
3 NO2	0.000	mg/L
4 Br	0.000	mg/L
5 NO3	9.00	mg/L
6 PO4	0.000	mg/L
7 SO4	30.8	mg/L




Cation		
Factor		
1 Li	0.000	mg/L
2 Na	9.55	mg/L
3 NH4-N	0.000	mg/L
4 K	2.39	mg/L
5 Mg	4.49	mg/L
6 Ca	23.0	mg/L




### ■ Liquid fertilizer

Anion		
Factor		
1 F	1000	mg/L
2 Cl	16.0	mg/L
3 NO2-N	0.000	mg/L
4 Br	0.000	mg/L
5 NO3-N	1.20	g/L
6 PO4-P	4.44	g/L
7 SO4-S	79.0	mg/L




Cation		
Factor		
1 Li	10000	mg/L
2 Na	2.17	g/L
3 NH4-N	33.1	g/L
4 K	50.3	g/L
5 Mg	0.000	mg/L
6 Ca	0.000	mg/L
GH	0.000	mg/L




### ■ Rain water

Anion		
Factor		
1 F	0.000	mg/L
2 Cl	0.480	mg/L
3 NO2	0.018	mg/L
4 Br	0.000	mg/L
5 NO3	3.11	mg/L
6 PO4	0.000	mg/L
7 SO4	2.17	mg/L




Cation		
Factor		
1 Li	0.000	mg/L
2 Na	0.330	mg/L
3 NH4	0.994	mg/L
4 K	0.336	mg/L
5 Mg	0.052	mg/L
6 Ca	0.935	mg/L
GH	2.55	mg/L




### ■ Wastewater

Anion		
Factor		
1 F	5	mg/L
2 Cl	0.085	mg/L
3 NO2-N	0.000	mg/L
4 Br	0.000	mg/L
5 NO3-N	16.3	mg/L
6 PO4	1.86	mg/L
7 SO4	31.3	mg/L




Cation		
Factor		
1 Li	0.000	mg/L
2 Na	45.6	mg/L
3 NH4-N	0.078	mg/L
4 K	10.8	mg/L
5 Mg	5.79	mg/L
6 Ca	25.5	mg/L
GH	87.4	mg/L




## ●Examples of three types of hardness printing

### ■ Mineral water Domestic

Cation		
Factor		
1 Li	0.000	mg/L
2 Na	4.42	mg/L
3 NH4	0.000	mg/L
4 K	0.531	mg/L
5 Mg-hd	8.33	mg/L
6 Ca-hd	17.0	mg/L
GH	25.4	mg/L




Cation		
Factor		
1 Li	0.176	mg/L
2 Na	78.4	mg/L
3 NH4	0.000	mg/L
4 K	7.54	mg/L
5 Mg-hd	120	mg/L
6 Ca-hd	408	mg/L
GH	528	mg/L



### ■ Deep water


Cation		
Factor		
1 Li	0.000	mg/L
2 Na	18.1	mg/L
3 NH4	0.000	mg/L
4 K	15.2	mg/L
5 Mg-hd	208	mg/L
6 Ca-hd	48.9	mg/L
GH	255	mg/L




## ●Example of printing with total hardness only

### ■ Beer

Anion		
Factor		
1 Li	0.480	mg/L
2 Na	23.7	mg/L
3 NH4	14.6	mg/L
4 K	259	mg/L
5 Mg	88.5	mg/L
6 Ca	39.6	mg/L
GH	381	mg/L




Anion		
Factor		
1 Li	0.000	mg/L
2 Na	24.8	mg/L
3 NH4	0.820	mg/L
4 K	156	mg/L
5 Mg	36.7	mg/L
6 Ca	34.8	mg/L
GH	238	mg/L




## ●Others

### ■ Spring water

Cation		
Factor		
1 F	0.067	mg/L
2 Cl	9.34	mg/L
3 NO2	0.000	mg/L
4 Br	0.000	mg/L
5 NO3	5.24	mg/L
6 PO4	0.000	mg/L
7 SO4	69.6	mg/L



Anion		
Factor		
1 Li	0.000	mg/L
2 Na	7.89	mg/L
3 NH4	0.000	mg/L
4 K	1.05	mg/L
5 Mg	10.2	mg/L
6 Ca	51.0	mg/L





# Peripheral equipment

Reduced labor by automatic measurement of multiple samples

## Autosampler ICA-700AS

### ■ Features

#### ● Simply connect with a dedicated cable\*1

When it is connected by a dedicated cable, the sample in the vial of No.1 is judged as a standard solution, and automatic calibration is performed with IA-300.

#### ● Automatic measurement of 100 continuous samples

The vial rack is equipped with two 50 samples as standard.\*2  
(One sample is used as a standard solution.)

Clean the aspiration needle after each sample measurement.\*3

#### ● Interaction with IA-300 Timer Function

IA-300 can be automatically started and stopped



Connecting to an IA-300



### ■ Specifications

Model name	ICA-700AS
Display	Backlit LCD
Sample container	1.5mL dedicated container Vial(100 pieces): 136C408 Vial cap(100 pieces): 136C409
Injection volume of sample	1 to 1000 $\mu$ L (1 $\mu$ L step)
Number of samples	Up to 100 samples (50 samples x 2)
Sample injection type	Loop mode
Material of wet part	PEEK, ETFE, $\beta$ titanium
Operating temperature range	4 to 35°C
Cooler Unit (Optional) *4	Presence/absence
Power supply	AC100 to 240V 50/60Hz
Power consumption	Max. 20VA
External dimensions	263(W) × 220(H) × 416(D)mm
Weight	11.8kg



With cooler unit (optional)

\*1 Please specify IA-300 as the connected model when ordering. If the connection model is not specified, the connection cable is not included.

The code for purchasing the connection cable separately is 7678110K.

\*2 Vial racks can be purchased additionally for sample preparation.

\*3 Use ultrapure water as cleaning water.

\*4 Cooler unit is sold separately. The code for purchasing separately is 133A047.

# Peripheral equipment

## For data management using a PC

### Data-gathering software **GP-LOG**

The data can be imported to a PC in a textual format by RS-232C.

The software can be downloaded free of charge after you register.

RS-232C cable is sold separately. (Code Number: 0GC00002 )

### ■ Convenient function

**The results can be stored after the measurement is completed.**

Press the "Start" button on GP-LOG to wait while receiving, and press the "Stop" button after reading is completed to add the time-in period and send the data to the PC.

	A	B	C	D	E	AE	AF	AG	AH	AI	AJ	AK	AL	AM	AN	AO	AP	AQ	AR
1	2017/5/9	13:02:19	DM	1	2017/4/12 14:32	2.01	mg/L	Mg	1	7.71	mg/L	Ca	1	20.5	mg/L	GH	1	83.1	mg/L

**You can load data (up to 100 data) stored in IA-300 memory.**

Starting No. and ending No. can be set from the "DATA" key on IA-300 and can be imported into the PC. (For 1-data, set the starting No. and ending No. to the same.)

	A	B	C	D	E	AE	AF	AG	AH	AI	AJ	AK	AL	AM	AN	AO	AP	AQ	AR
1	2017/5/9	13:03:07	DM	1	2017/4/12 14:32	2.01	mg/L	Mg	1	7.71	mg/L	Ca	1	20.5	mg/L	GH	1	83.1	mg/L
2	2017/5/9	13:03:07	DM	2	2017/4/12 15:35	0.949	mg/L	Mg	1	1.05	mg/L	Ca	1	4.56	mg/L	GH	1	15.7	mg/L
3	2017/5/9	13:03:07	DM	3	2017/4/12 16:07	0.047	mg/L	Mg	1	0.917	mg/L	Ca	1	8.43	mg/L	GH	1	24.8	mg/L
4	2017/5/9	13:03:08	DM	4	2017/4/12 16:28	0.138	mg/L	Mg	1	0.909	mg/L	Ca	1	1.24	mg/L	GH	1	6.83	mg/L
5	2017/5/9	13:03:08	DM	5	2017/4/13 9:18	5.17	mg/L	Mg	1	2.32	mg/L	Ca	1	5.11	mg/L	GH	1	22.3	mg/L
6	2017/5/9	13:03:08	DM	6	2017/4/13 9:38	0	mg/L	Mg	1	0.486	mg/L	Ca	1	0.068	mg/L	GH	1	2.17	mg/L
7	2017/5/9	13:03:08	DM	7	2017/4/13 10:06	3.96	mg/L	Mg	1	1.31	mg/L	Ca	1	0.143	mg/L	GH	1	5.76	mg/L
8	2017/5/9	13:03:08	DM	8	2017/4/13 10:28	2.54	mg/L	Mg	1	0.002	mg/L	Ca	1	0.581	mg/L	GH	1	1.46	mg/L
9	2017/5/9	13:03:09	DM	9	2017/4/13 10:57	2.53	mg/L	Mg	1	0.448	mg/L	Ca	1	0.946	mg/L	GH	1	4.21	mg/L
10	2017/5/9	13:03:09	DM	10	2017/4/13 11:54	1.14	mg/L	Mg	1	0.636	mg/L	Ca	1	0.302	mg/L	GH	1	3.37	mg/L

**Real-time data acquisition is also possible.**

Press the "Start" button on GP-LOG to wait while receiving, and after completing measurement, press the "Stop" button for each measurement to add the take-in period and send the data to the PC.


	A	B	C	D	E	F
1	2017/5/9	13:04:28	CD		7	
2	2017/5/9	13:04:28	100 0			
3	2017/5/9	13:04:28	101 0000.000	02524.638	000.0	08.5 39 39.7
4	2017/5/9	13:04:28	101 0000.100	02524.638	000.0	08.5 39 39.7
5	2017/5/9	13:04:28	101 0000.200	02524.638	000.0	08.5 39 39.7
6	2017/5/9	13:04:28	101 0000.300	02524.638	000.0	08.5 39 39.7
7	2017/5/9	13:04:28	101 0000.400	02524.638	000.0	08.5 39 39.7
8	2017/5/9	13:04:28	101 0000.500	02524.638	000.0	08.5 39 39.7
9	2017/5/9	13:04:28	101 0000.600	02524.638	000.0	08.5 39 39.7
10	2017/5/9	13:04:29	101 0000.700	-2480.477	000.0	08.5 39 39.7
11	2017/5/9	13:04:29	101 0000.800	-2480.477	000.0	08.5 39 39.7
12	2017/5/9	13:04:29	101 0000.900	-2480.477	000.0	08.5 39 39.7
13	2017/5/9	13:04:29	101 0001.000	-2480.477	000.0	08.5 39 39.7
14	2017/5/9	13:04:29	101 0001.100	-2480.477	000.0	08.5 39 39.7
15	2017/5/9	13:04:29	101 0001.200	-2480.477	000.0	08.5 39 39.7

10790	2017/5/9	13:22:26	101 1078.700	-0019.569	000.0	08.5 39 39.7
10791	2017/5/9	13:22:26	101 1078.800	-0019.569	000.0	08.5 39 39.7
10792	2017/5/9	13:22:26	101 1078.900	-0019.552	000.0	08.5 39 39.7
10793	2017/5/9	13:22:27	101 1079.000	-0019.542	000.0	08.5 39 39.7
10794	2017/5/9	13:22:27	101 1079.100	-0019.542	000.0	08.5 39 39.7
10795	2017/5/9	13:22:27	101 1079.200	-0019.544	000.0	08.5 39 39.7
10796	2017/5/9	13:22:27	101 1079.300	-0019.544	000.0	08.5 39 39.7
10797	2017/5/9	13:22:27	101 1079.400	-0019.556	000.0	08.5 39 39.7
10798	2017/5/9	13:22:27	101 1079.500	-0019.572	000.0	08.5 39 39.7
10799	2017/5/9	13:22:27	101 1079.600	-0019.572	000.0	08.5 39 39.7
10800	2017/5/9	13:22:27	101 1079.700	-0019.589	000.0	08.5 39 39.7
10801	2017/5/9	13:22:27	101 1079.800	-0019.589	000.0	08.5 39 39.7
10802	2017/5/9	13:22:27	101 1079.900	-0019.605	000.0	08.5 39 39.7
10803	2017/5/9	13:22:28	103			
10804	2017/5/9	13:22:28	102 0			




# Setting GP-LOG of the data-collection software

## ■ Setting of the equipment used



Select etc. from the list of devices to use.

## ■ Setting the conditions for compilation



Select receive only.

## Specifications

Item	Contents
Model name	IA-300
Measurement method	Ion chromatography
Measurement items	Anion (non-suppressor type) PO <sub>4</sub> , F, Cl, NO <sub>2</sub> , Br, NO <sub>3</sub> , SO <sub>4</sub> Anion (suppressor type) F, Cl, NO <sub>2</sub> , Br, NO <sub>3</sub> , PO <sub>4</sub> , SO <sub>4</sub> Cation 1 and divalent simultaneous Li, Na, NH <sub>4</sub> , K, Mg, Ca
Repeatability	2% C.V. or less in the calibration solution
Sample injection	Manual sample injection and manual valve switching
Sample measurement	Loop-cut loop volume 20 $\mu$ L or 200 $\mu$ L
Measurement time	15 to 18 min/batch (according to measurement conditions) * Analysis is possible in 10 minutes only when using the old column PCI-302S column.
Calibration	One-point calibration using the specified calibration solution
Column oven	40 $\pm$ 4°C
Data Processing	Built in
Detection part	Method : Electrical conductivity detection Cell temperature control : 40 $\pm$ 4°C
Display	Graphic LCD
Printer	Built-in thermal printer
Operating temperature range	10 to 35°C However, no abrupt temperature change
Output	Analogues: 0 to 1V Digital: RS-232C
Power supply	AC100V 50/60Hz
Power consumption	Max. 250VA
Dimensions and weight	Approx. 190(W) $\times$ 469(H) $\times$ 530(D)mm and approx. 18kg

## Standard accessories


- 1mL disposable syringe
- Syringe needle
- Sample loops (20, 200  $\mu$ L) (1 each)
- Syringe set for air bleeding
- Spanners (6 x 8, 8 x 10) (1 each)
- Hexagon wrenches (1.5mm, 2.5mm, 3mm) (1 each)
- Plunger seal replacement jig
- Printer paper (Volume 2)
- AC cable
- 2P transmitter adapter
- Ground wire
- Ring instructions manual

\* Column, chemical suppressor, calibration solution, lysis solution, removal solution and tank introduction pipe are sold separately.

The required parts differ depending on the measurement mode, so please select from "Required parts for each measurement mode" on page 10.

## Measurement range

Mode/Measurement Ion Type		When 20 $\mu$ L loop is used	When using a 200 $\mu$ L loop
PCI-322 / 1, 2 divalent cation determination	Li	0.050 to 10.00mg/L	0.005 to 1.00mg/L
	Na, Mg, NH <sub>4</sub> (NH <sub>4</sub> -N)	0.250 to 50.0mg/L (0.194 to 38.8mg/L)	0.025 to 5.00mg/L (0.019 to 3.88mg/L)
	K, Ca	0.500 to 100mg/L	0.050 to 10.0mg/L
PCI-201S / Anion measurement (non-suppressor)	F, Cl, Br	1.00 to 100mg/L	0.100 to 10.0mg/L
	NO <sub>2</sub> (NO <sub>2</sub> -N)	1.00 to 100mg/L (0.305 to 30.5mg/L)	0.100 to 10.0mg/L (0.031 to 3.05mg/L)
	NO <sub>3</sub> (NO <sub>3</sub> -N)	1.00 to 100mg/L (0.226 to 22.6mg/L)	0.100 to 10.0mg/L (0.023 to 2.26mg/L)
	SO <sub>4</sub> (SO <sub>4</sub> -S)	2.00 to 200mg/L (0.668 to 66.8mg/L)	0.200 to 20.0mg/L (0.067 to 6.68mg/L)
	PO <sub>4</sub> (PO <sub>4</sub> -P)	5.00 to 200mg/L (1.63 to 65.2mg/L)	0.500 to 20.0mg/L (0.163 to 6.52mg/L)
PCI-211 / Anion measurement (non-suppressor)	F, Cl, Br	0.500 to 50.0mg/L	0.050 to 5.00mg/L
	NO <sub>2</sub> (NO <sub>2</sub> -N)	0.500 to 50.0mg/L (0.152 to 15.2mg/L)	0.050 to 5.00mg/L (0.015 to 1.52mg/L)
	NO <sub>3</sub> (NO <sub>3</sub> -N)	0.500 to 50.0mg/L (0.113 to 11.3mg/L)	0.050 to 5.00mg/L (0.011 to 1.13mg/L)
	SO <sub>4</sub> (SO <sub>4</sub> -S)	1.00 to 100mg/L (0.334 to 33.4mg/L)	0.100 to 10.0mg/L (0.033 to 3.34mg/L)
	PO <sub>4</sub> (PO <sub>4</sub> -P)	2.50 to 100mg/L (0.815 to 32.6mg/L)	0.250 to 10.0mg/L (0.082 to 3.26mg/L)
PCI-205 / Anion measurement (suppressor)	F, Cl, Br	0.050 to 50.0mg/L	
	NO <sub>2</sub> (NO <sub>2</sub> -N)	0.050 to 50.0mg/L (0.015 to 15.2mg/L)	
	NO <sub>3</sub> (NO <sub>3</sub> -N)	0.050 to 50.0mg/L (0.011 to 11.3mg/L)	
	SO <sub>4</sub> (SO <sub>4</sub> -S)	0.100 to 100mg/L (0.033 to 33.4mg/L)	
	PO <sub>4</sub> (PO <sub>4</sub> -P)	0.250 to 100mg/L (0.082 to 32.6mg/L)	





## Required parts for each measurement mode (sold separately)


Measurement mode	Parts		Necessary items	Quantity, etc.
	Name	Model name		
Anion measurement PCI-205 mode Suppressor	Anion exchange column	PCI-205	○	1 piece
	Guard column	PCI-205G	○	1 piece
	Chemical suppressor	6813690K	○	1 piece
	Calibration solution	For 20 $\mu$ L loop IA-AS1	○	100mL 1 bottle
	Eluent	1L 143H063	○	1 bottle
	Removal solution	1L 143H071	○	1 bottle
Anion measurement PCI-211 mode Non-suppressor	Tank inlet pipe	6547830K	○	With 2L tank
	Anion exchange column	PCI-211	○	1 piece
	Guard column	PCI-211G	○	1 piece
	Calibration solution	For 20 $\mu$ L loop IA-AS1	△*1	100mL 1 bottle
		For 200 $\mu$ L loop IA-AS2	△*1	100mL 1 bottle
	Eluent	2L 6547760K	○	1 bottle
		5L 6547770K	△*2	1 bottle
Anion measurement PCI-201S mode Non-suppressor	Tank inlet pipe	6547830K	○	With 2L tank
	Anion column	PCI-201S	○	1 piece
	Guard column	PCI-201SG	○	1 piece
	Calibration solution	For 20 $\mu$ L loop IA-AS1	△*1	100mL 1 bottle
		For 200 $\mu$ L loop IA-AS2	△*1	100mL 1 bottle
	Eluent	2L IA-AE-12	○	1 bottle
		5L IA-AE-15	△*2	1 bottle
Cation measurement PCI-322 mode Simultaneous determination of 1 and 2 valences	Tank introduction pipe	6547830K	○	With 2L tank
	Cation column	PCI-322	○	1 piece
	Guard column	PCI-322SG	○	1 piece
	Calibration solution	For 20 $\mu$ L loop IA-CS1	△*1	100mL 1 bottle
		For 200 $\mu$ L loop IA-CS2	△*1	100mL 1 bottle
	Eluent	2L 143H061	○	1 bottle
		5L 143H062	△*2	1 bottle
Cation measurement PCI-322 mode Simultaneous determination of 1 and 2 valences	Tank inlet pipe	6547830K	○	With 2L tank

\*1 : Requires either a 20  $\mu$ L loop or a 200  $\mu$ L loop. \*2 : 5L is used for refilling.

## Other common parts


Parts		Quantity	Notes
Name	Model/Code No.		
Printer paper	PAP-HCS	5 Volume	Thermal recording paper
RS-232C connecting cable	0GC00002	1 piece	A commercially available USB serial converter is required to connect to USB.
IA-300 ICA-700AS connection cable	7678110K	1 piece	For connecting autosampler ICA-700AS

## Dimensions (Unit :mm)



**We provide high-level solutions to analytical needs in a wide range of fields, including the environment, effluent, agriculture, food, water, quality control, and education.**

## Application areas



### Agriculture

Soil, hydroponic solution, fertilizer, livestock feed, vegetables, tea leaves, etc.




### Wastewater

Plant effluent, semiconductor plant effluent, waste water from waste treatment plants, sewage, etc.



### Tap water

Raw water, purified water, etc.




### School

Support for research, education, etc.




### Analysis

For request analysis work




### Environmental

Rainwater, rivers, lakes, groundwater, soil, rocks, seawater, etc.



### Food products

Soft drinks, juices, natural waters, liquors, foods, food additives, etc.



### Manufacturing

Raw materials, washing water, products, plating solution, treatment solution, etc.



### Others

Various research, quality control, etc.



**DKK-TOA CORPORATION**

Overseas Sales Division:  
DKK-TOA Corporation  
29-10, 1-Chome, Takadanobaba, Shinjuku-ku,  
Tokyo 169-8648 Japan  
Tel : +81-3-3202-0225 Fax : +81-3-3202-5685  
E-mail : intsales@dkktoa.com



**CAUTION**

Please read the operation manual carefully before using products.