

SAT-500 SALT ANALYZER

TORDER

換算值

AT IT'S THE

1.008

1.008

SALT AMALYZER SAT-500

STOP

TOADK

Quick, Easy, Extremely accurate

Expert salinity measurement for quality control on food production lines

The Model SAT-500 provides easy measurement of salinity (chloride concentration) by simply dispensing a small amount of the sample solution using a dispenser. In addition, the measurement is quick, with the time from the start of measurement to data output within about 25 seconds.

It is equipped with a variety of functions suitable for use in a wide range of applications, such as quality control in food processing, as well as labor saving, time reduction, and high reliability in HACCP.

DKK-TOA CORPORATION

SAT-500 salt analyzer

Measurement Procedure



Features

Simple operation

The basic measurement operation is to inject a sample (solution) into the electrolyte (dedicated cell) using the dispenser.

Quick measurement

Measurement and data output are completed in about 25 seconds for 1% standard solutions.

Extremely accurate measurement

Using the 20 $\,\mu\text{L}$ dispenser, a measurement with less than 0.5% in C.V. is possible for 1% NaCl standard solution.

Same electrolyte (dedicated beaker) can be used continuously

The same electrolyte can be used for continuous measurement; more than 100 measurements for 1% NaCl (when using a 20 μ L dispenser) or about 50 measurements for soy sauce (when using a 2 μ L dispenser).

Save running cost and minimize routine set-up work

The same electrolyte can be used for continuous measurement. The same solution can be used for more than 100 measurements of 1% standard solution or 50 measurements of

high concentration samples.

This reduces running costs and minimizes routine set-up work.



DKK-TOA CORPORATION

Overseas Sales Division:

DKK-TOA Corporation

29-10, 1-Chome, Takadanobaba, Shinjuku-ku, Tokyo 169-8648 Japan Tel: +81-3-3202-0225 Fax: +81-3-3202-5685

Principle of Measurement

The Model SAT-500 measures the chloride content to obtain salinity based on coulometric titration method (end point is detected by potential difference). When an electric current passes between the silver wire(SAT-AG) and the silver electrode(AG-111), white deposits are formed from the dissolved silver ions and the chloride present in the sample: $Ag^+ + CI^- \rightarrow AgCI \downarrow$

Therefore, the chloride content in the sample can be determined from the electric charge (current x time) required until all of the chloride in the sample has



been depleted by electrolysis. Since the Model SAT-500 uses a highly stable constant current for electrolysis, the electric charge is proportional to time, enabling it to obtain the chloride content by measuring the time needed for electrolysis. In addition, the end point is detected by sensitive calculation of the time when chloride ions disappear and surplus silver ions start to appear using the silver electrode (AG-111) and reference electrode (ELR-003).

Specifications

-			
Measuring Method	Coulometric titration method (end point is detected by potential difference).		
Display	Graphic LCD panel (320 x 240 pixels) with back light. Selectable English/ Japanese language menus.		
Sample Injection Volume:	2, 20, 200μL or OTHER (1 ~ 9999μL)		
Measuring Ranges (displayed by %)	2µLdispenser, 0.00 ~ 49.99% lower detection limit 2.0% 20µL dispenser, 0.000 ~ 4.999% lower detection limit 0.2% 20µL dispenser, 0.0000 ~ 0.4999% lower detection limit 0.02% 0THER disp. vol., 0.0000 ~ 99.99% (depends on calibration solution concentration)		
Concentration Conversion	NaCI/CI [®] switching		
Concentration Units	% (W/V), % (W/W), mg/L, mol/L		
Repeatability	0.5% or less in C.V. value at 1% of NaCl standard solution with 20 μL dispenser		
Calibration	Automatic single point calibration using NaCl standard solution at the concentration for the injection volume		
Measurement Start	Automatic start on detection of potential change or manual start by pressing START key.		
Printer	Built-in thermal printer		
Printer Data Memory	300 data points (with backup)		
Output	RS-232C (for PC or external printer)		
Ambient Temperature	$10 \sim 35^{\circ}$ C (avoid sudden temperature changes)		
Power Source	100 ~ 240VAC, 50/60Hz		
Power Consumption	Approx. 20VA (maximum)		
Dimensions, Weight	Main; 190(W) × 192(H) × 385(D)mm, Approx. 3.3 kg		
	Meas. Section; 110(W) × 370(H) × 220(D) mm, Approx. 1.1 kg		

Standard accessories

Product name	Model	Qty	Product name	Model	Qty
Stirrer ASSY	6877950K	1	Beaker for SAT 20mL,	SAT-B20	2
Electrode holder ASSY	6877940K	1	SAT-B2U		_
Silver electrode	AG-111	2	Stirrer (two), SAT-SB	SAT-SB	2
Silver wire for SAT	RATAC	1	Printer paper (2 rolls) %1	—	1
Silver wire für SAT	SAT-AG	1	Power code	1180025	1
Reference electrode	ELR-003	1		1100020	
Electrolyte for Salt Analyzer	143A373	1	2P Conversion adapter	102A991	1
SAT-1F1 500mL			Earth lead	6545820K	1
Addition Agent for	143A374	1	Instruction manual	145A474	1
Salt Analyzer SAT-TGT SUML			Simple instruction manual	145A475	1
Abrasive Solution for	143A376 143A375	1	Mioro disponsor #2	(specify conc.)	1
Salt Analyzer SAT-TZT 50mL			MICIO disperiser %2		
Reference Electrode Inner Solution for Salt Analyzer SAT-1E1 100mL			Standard Solution for Salt Analyzer %3		1

*1) When ordering additional printer paper, please order 5 rolls (PAP-HCS).
 *2) Only one set of micro-dispensers is included as part of standard accessories. Please select from 2 μL (OSZ00018), 20 μL(OSZ00021), 200 μL(OSZ00024).

#3) When ordering, please select from 10% NaCl Standard solution SAT-1A1 (143A371) or 1% NaCl Standard solution SAT-1A2 (143A372).

Options

 Data compilation software G-LOG2 •RS-232C cable (2m) 0GC00002

Micro Dispenser Components

•For 2µL : Glass tube (100 pcs) 0SZ00019, Parts kit for replacement 0SZ00020
 For 20μ
 Glass tube (100 pcs) 05Z00025, Parts kit for replacement 05Z0023

 For 20μL
 Glass tube (100 pcs) 05Z00025, Parts kit for replacement 05Z0026

* Breakdown of parts

3 stailess steel piston rod, 2 O-rings, 1 knob ring, 1 coil spring, 1 barrel tube, 3 Teflon chips (for 200 μ L only)

CAUTION

