SPECIFICATION SHEET

pH Electrode for Bioprocess DO Electrode for Bioprocess



ELP-045/046/047/048 ELD-011/012/013/015

The pH electrode and dissolved oxygen (DO) electrode for bioprocessing (fermentation) must be sterilized repeatedly with each container, and therefore must be constructed to withstand high pressure at high temperature.

Our bioprocess electrode exhibits excellent long-term performance even under such severe conditions of use.

Since both electrodes are lead-free, they are environmentally friendly.

Features of pH electrodes

- ○The glass electrode is safe to handle because the glass electrode uses a Strong film with mechanical strength about 10 times higher than that of our conventional glass film.
- ○The reference electrode is a KCl-free gel type that is easy to maintain. In addition, because of the double junction, the potential is stabilized and the silver ion trap (chelating resin) is used. The junction is hard to clog.
- OThe cooking sterilization temperature corresponds to 135°C.
- ○The mini-jar insertion type has four lengths. Stainless steel holders and welded sockets are also available for mounting on the side of the tank.

Features of Dissolved Oxygen (DO) electrodes

- ○The membrane of a polarographic electrode, which was previously a single-layer FEP membrane have realized with a SUS mesh use of sandwiched high-strength laminated films for improved durability and longer life.
- ○Since the internal electrode (working electrode and counter electrode) was made into a cartridge for the first time (Japanese Patent No. 4405653), only the electrodes can be replaced, reducing running costs.
- ○Since the electrode body is made of SUS316L, the electrode body can be handled safely.
- $\odot\ensuremath{\mathsf{There}}$ are no pressure balancing holes, so contamination does not occur.
- ○The mini-jar insertion type has three lengths. In addition, tank side mounting type and welded socket type are provided.



pH/DO electrode for the mini-jar



pH detector for the tank



Dissolved Oxygen (DO) detector for the tank



System configuration

pH Electrode for Bioprocess(Fermentation)

• For mini-jar ELP-045/046/047/048

Standard	Specifications
----------	----------------

Measurement method	method : Glass electrode method					
measurement range	measurement pH 0 1	measurement pH 0 to 12)				
Use Temperature	: Measured at 0 to 50°	ъС				
	Sterilized at 135°C o	or less				
Use Pressure	: Measured at 37°C 0.	5MPa or les	ss			
	Sterilized at 135°C ().3MPa or le	ess			
Internal electrode	: Silver-silver chloride	è				
Inner Solution for	: 3.3mol/L KCl gels					
Reference Electrode	<u>!</u>					
Material of wetted parts	s : glass, ceramics	glass, ceramics				
Lead wire	Leadless type, Connector connection					
Outer diameter of	f :φ12					
electrode glass part						
Mounting method	Mounting method : Insertion type					
Dimensions under	: As shown in the Model L (mm					
thread	right figure	ELP-045	120			
		ELP-046	220			
		ELP-047	320			

ELP-048

420





Mounting example for mini-jar

• Fixed holder for tank HEW-101

Standard Specifications

Combined Electrode	:	ELP-045
Material of wetted parts	: ;	SUS316L
Mounting flange	: (φ38
outer diameter		
Insertion part outer	: ($\varphi 24.5$
diameter		
O ring for seal	:	P21 FKM
Size below flange	:	116mm





• Steaming test data



Change in electromotive force of the $\ensuremath{\mathsf{pH7}}$





* Steaming conditions: Steaming is performed 50 times at 121°C for 30 minutes.

• The special cable ELW-046/047/048/049

Standard Specifications



Mode	L L (m)
ELW-04	16 1
ELW-04	17 3
ELW-04	18 6
ELW-04	19 10

DO Electrode for bio-process (fermentations)

• For mini-jars ELD-011/012/013





Schematic cross-sectional view of the laminated film

• Standard accessories

Model	Name	Amount
	Exchange membrane for $\varphi 12$ (with O-ring)	5
	Electrolyte 50mL	1
	Syringe (1mL)	1
ELD-015	Blind cap	1
	O-ring (P10A)	1
	Exchange membrane for $\varphi 25$ (with O-ring)	5
	Electrolyte 50mL	1
ELD-015	Syringe (2.5mL)	1
	Blind cap	1
	O-ring (P10A)	1

Supplements

Model	Name	code No.	Amount
	Exchange membrane for $\varphi 12$ (with O-ring)	0CA00001	5
	Electrolyte 50mL	143A461	1
ELD-012 ELD-013	Electrode cartridges (for φ12) ELD-027	ELD027	1
	Sodium Sulfite	143A030	1
	Exchange membrane for φ25 (with O-ring)	0CA00002	5
ELD-015	Electrolyte 50mL	143A461	1
	Electrode cartridges (for φ25) ELD-028	ELD028	1
	Sodium Sulfite	143A030	1

• Steaming test data



*Steaming condition: Steaming at 121°C for 30 minutes is performed 20 times without replacing the diaphragm and Electrolyte.

The special cable ELW-074

Standard Specifications



Mounting Equipment

Equipment to attach to the tank for the pH detector HEW-101 type (built-in pH Electrode ELP-045 type) and DO (DO) detector ELD-015 type for mounting to the side of the tank.

Welded sockets, cap nuts, blind plugs, etc. are provided.

• Welded socket (for tank) Code No.: X0134900





• Weld socket cap nut Code No.: X0110100



G1 1/4

Blind plug (for tank) Code No.: BP-0149

This plug is inserted into the socket of the tank when the detector is removed.



Calibration equipment

Calibration kits for calibrating pH meters and dissolved oxygen (DO) meters for bioprocesses, and portable pH meters and portable dissolved oxygen (DO) meters for cross-checking are available.

pH calibration kit

The calibration kit for the pH meter consists of a pH7 reference solution for two-point calibration and a pH4 or pH9 reference solution and beaker. Select pH4 standard solution kit for standard measurement and pH9 standard solution kit for alkaline measurement.

Applications	Products Name	Code No.
pH Standard	pH Calibration 4-7 Kit (AN)	$6531680 { m K}$
pH alkali	pH calibration 7-9 kit (BN)	6531690 K

Kit breakdown table

○ pH Calibration 4-7 Kit (AN)	Code No. Q	uantity
Polyethylene beaker (500mL)	(136C035)	3
pH4 standard solution (500mL)	(143F191)	1
pH7 standard solution (500mL)	(143F192)	1
Powder reagent for pH4 reference solution	(143F060)	1
(5 bags/set for 500mL)		
Powder reagent for pH7 reference solution	(143F061)	1
(5 bags/set for 500mL)		

\odot pH calibration 7-9 kit (BN)	Code No. Quantity	
Polyethylene beaker (500mL)	(136C035)	3
pH7 standard solution (500mL)	(143F192)	1
pH9 standard solution (500mL)	(143F193)	1
Powder reagent for pH7 reference solution	(143F061)	1
(5 bags/set for 500mL)		
Powder reagent for pH9 reference solution	(143F062)	1
(5 bags/set for 500mL)		

Dissolved Oxygen (DO) Calibration Kit

There are two calibration kits for the DO meter: one is for span calibration with air saturated water (A) and the other is for span calibration in air (B).

Product name	Applications	Code No.
DO calibration kit (A)	For span calibration with air saturated water	6151840K
DO calibration kit (B)	For span calibration in air	6151850 K

[Contents of DO calibration kit (A)]

Item No.	Product name	Code No.	Quantity	Notes	Applications
1	Special grade anhydrous sodium sulfite	143A002	1	500 g bottle	For preparing zero water
2	Beaker	136C035	3	500mL made of PE	For calibration and cleaning
3	Air pump	125B248	1	Dry cell type	For air saturated water preparation
4	Dry cell	106A003	2	Single 1 1.5V	For air pump

[Contents of DO calibration kit (B)]

Item No.	Product name	Code No.	Quantity	Notes	Applications
1	Special grade anhydrous sodium sulfite	143A002	1	Contains 500g	For preparing zero water
2	Beaker	136C035	3	500mL made of PE	For calibration and cleaning

Por table pH meter



• Portable dissolved oxygen (DO) meter

For DO cross check	
DO-31P	
A CONTRACT OF A	
Power-saving design, rechargeable batteries also availab	le
• Waterproof construction (IP67:1m, immersion allowed	
for 30minutes)	

- Data memory function (1000 data)
- Measurement range 0 to 20.00 mg/L,0 to 200%, 0 to 50°C

Combined pH transmitter



- temperature) function •One-touch automatic calibration, various self-
- diagnosis functions
- ●Meas. range:pH -1 to 15

TON DIKK

Overseas Sales Division: DKK-TOA Corporation 29-10, 1-Chome, Takadanobaba, Shinjuku-ku, Tokyo 169-8648 Japan Tel : +81-3-3202-0225 Fax : +81-3-3202-5685 E-mail : intsales@dkktoa.com



Combined DO transmitter





•One-touch automatic calibration, various selfdiagnosis functions

before using producuts.

Please read the operation manual carefully

• Meas. range: Minimum : 0 to 1 mg/L Maximum : 0 to 50 mg/L

https://www.toadkk.com/english/

CAUTION